

SURCAM

Dental Pro.

Cleaning and Sterilization instructions

Product features:

A. Description:

1. Closed-tray Transfers, Healing Caps made of titanium Ti6Al4V-Eli for multiple use.
2. Esthetic Abutment and Multi-Units is delivered non-sterile made of titanium Ti6Al4V-Eli for single use and need to be sterilized prior to use.
3. Temporary abutments made of PEEK and need to be sterilized prior to use.
4. Open tray transfers and accessories are made of SAE 303 stainless steel and need to be sterilized prior to use.

B. Packaging features:

Items are packaged in bags suited to preserve functional device characteristics and cleanliness.

C. Warnings:

All instruments must be cleaned (ultrasound equipment is recommended) and sterilized before use as per regulation UNI EN 556-1:2002.

Cleaning and sterilization must only be performed with specific materials, specifically detergent containing the following substances should be avoided:

- Oxalic acid.
- High concentrate chlorine.

Insert the instruments in the disinfectant solution immediately after surgery to prevent encrustation of blood, secretions, etc.

Do not store wet or damp instruments.

Do not sterilize, clean or disinfect instruments made with different metals in the same cycle.

Steam sterilization in autoclave is specifically recommended.

D. Product identification:

SURCAM'S Products are identified by the item code and production lot printed on the packaging.

E. Specific liability:

SURCAM'S Oral implant products meet the essentials requisites set by directive 93/42/CEE and directive 2007/47/CE.



Any product use other than that specified is considered "improper" and releases the manufacturer from any liability.

Incorrect or insufficient maintenance may damage instruments over the short-term.

Specific instrument cleaning products are readily available, strictly follow their instructions.

Warning: SURCAM's device use is reserved to dental surgeons. For information, clarifications or suggestions please contact SURCAM.

Packaging symbol legend:

<p>SURCAM <i>Dental Pro.</i></p> <p>Northern Industrial Zone P.O.B 12084 Nahariya 2201202, Israel Technical support: 972-4-9523334 Sales: 972-4-9523511 www.surcamdental.com E-mail: sales@surcamdental.com</p>	→	Manufacturer's Adress
<p>Abutment - finale 3.25 Q.ta:01</p>	→	Description and quantity
<p>REF SMS-0040</p>	→	Article code
<p>LOT 14-000088</p>	→	Production lot
<p>CE 0483</p>	→	CE mark
<p></p>	→	See the manual
<p></p>	→	Not sterile

Devices: All the surgical instruments manufactured by SURCAM included fixed or moving parts.

Instructions

Use way	Devices should be cleaned and sterilized before use.
Cleanliness preparation	All the devices with moving parts must be disassembled. It is better and recommended to make the preliminary cleanliness and the decontamination by Ultrasonic machine.
Cleaning by Ultrasonic machine	Immerse the instruments in an Ultrasonic bath filled with NDP detergent for 20 min. heated to 60°C temperature.
Manual cleaning	Wash the instruments under running tap water, removing detergent residues.
Sterilization	Using autoclave for 15 minutes for minimum of 132° C.
Storage	In a clean and dry place.
Manufactured contact	Very important is to know the phone number of the manufacture to solve possible problems.